



A TimeSens newsletter is born!

Dear TimeSens user,

You are reading the first issue of the TimeSens newsletter, introducing a new way to keep you informed about TimeSens evolution.

This first issue of the newsletter aims at introducing a new TimeSens version which has just been released.

This new version (0.95.0630) improves design, ergonomics and performances of TimeSens. It also includes several bug fixes, function enhancements and new features.

You can still use this version for free, until the commercial version is launched this autumn. The next newsletter will tell you everything about this commercial version.

Users having installed TimeSens before April 2014 must uninstall it, and then reinstall **TimeSens for Panel Leader** from www.timesens.com to access the new version. For the other users, this new version should install itself automatically.

Enjoy designing your innovative sensory sessions and analyzing the data collected with TimeSens.

As ever, your feedbacks on this new version are most welcome!


The TimeSens team




TimeSens Newsletter #1

What has changed?


- No more button bar, but contextual menus.
- Protocol, deployment, data and analysis tab contents displayed as tile views.
- New block view, experimental design view and upload view.
- More information displayed on products and descriptors views (related designs and labels).
- Fully customizable layout of scenario and output views.
- New simulation options in simulation view.




New protocol contextual menu & Tile view



New experimental design view



New scenario view



New simulation view

All these changes may have an impact on sessions developed in former versions, please double check your experimental designs before uploading those versions for deployment.



What's new?

- New “Recent files” menu.
- Custom fields for products.
- Multiple screen resolutions.
- Conditional visibility applied to single screens.
- New monitoring options (reset progress to specific screens, reset minimum time between connections, panelist browser information).
- Memorized analysis.
- New question table control.
- Conditional jump with questions.
- Import from XML session files.

What can you expect in future releases?

- Napping, Sorting, CATA protocols and analyses
- Improved tickets
- Tutorial (in French and English)
- Website update